


Patient Monitoring System Manual

- 1) Connect coordinator XBEE PRO S1 module to computer USB port.



- 2) Run Zigbee.exe program and select the suitable com port and baud rate. (Used COM3 and 9600 baud rate for this demo).



Form1

Com Port: COM3

Baud Rate: 9600

Host Address : http://127.0.0.1/

Connect

Disconnect

Received Data

TAG


BPM

Patient Details


Name :

Time :

- 3) Fill Host Address with system address and press Connect. (Demo using <http://p.site40.net>)



- 4) ON the router XBEE PRO s1



5) The data from router XBEE PRO s1 is shows at PMS system as below.



contact person:

Mohd Nur Fahmie

fahmie@semarakmotivasi.com

0134873831